# Cenefom

Cenefom Corp.

2F&3F., No.50-3, Keyan Rd., Hsinchu Science Park, Zhunan Township, Miaoli County 35053, Taiwan, R.O.C.

Tel: +886-37-580058

Fax: +886-37-580518

E-mail: [service@cenefom.com](mailto:service@cenefom.com) [www.cenefom.com](http://www.cenefom.com)

MedNet EC-REP GmbH

Borkstrasse 10, 48163 Münster, Germany

## Instructions for use

### ENT surgical sponge

(for the series NP-LS, NP-LL, NP-ML, NP-PL, NP-DL, NP-EL, NP-ES, NP-EB, NP-SL, EP-07, EP-09, EP-12 and EP-15)

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### Spugna Chirurgica ORL

(per le serie NP-LS, NP-LL, NP-ML, NP-PL, NP-DL, NP-EL, NP-ES, NP-EB, NP-SL, EP-07, EP-09, EP-12 e EP-15)

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### Esponja quirúrgica ORL

(para las series NP-LS, NP-LL, NP-ML, NP-PL, NP-DL, NP-EL, NP-ES, NP-EB, NP-SL, EP-07, EP-09, EP-12 y EP-15)

# English

## **1. Product Description:**

The ENT surgical sponge is provided sterile and for single use only. The device consists of PVA (polyvinyl alcohol) material which is capable of absorbing fluid several times of its own weight.

| Product Type   | Applicable surgery or treatment                |
|--|--|
| Standard Type<br>Large standard Type<br>Medium standard Type<br>Soft standard Type | Septoplasty, turbinectomy, rhinoplasty surgery |
| Epistaxis Type   | Nosebleed treatment                            |
| Sinus Type   | Sinus surgery                                  |
| Anatomical Type  | Nasal surgery to fill the nasal cavity         |
| Ear Pack   | Otic surgery                                   |
| Ear Wick   | Otitis externa related surgery                 |

## **2. Intended Use:**

The ENT surgical sponge is designed for placing inside the nasal cavities or ear canals as a tamponade or hemostasis material to control oozing or persistent bleeding.

## **3. Intended Therapeutic Use and Indications:**

The ENT surgical sponge is mainly used for hemostasis for ear canal (meatus), otitis chambers, bleeding within nasal cavity. Bleeding can be stopped effectively after insertion of ENT surgical sponge to provide a tamponade effect. It is intended to be used either after performing septoplasty, sinus surgeries, other nasal cavity surgeries and otology operations of the external auditory canal operations such as canalplasty or myringoplasty. An appropriate length and specification of polymer packing (or also containing a silicone vent tube) can be placed for nasal and otology surgeries at the bottom of the meatus or nasal cavity (single side or bilaterally).

## **4. Intended Patient Population:**

Male and female older than 7 years.

## **5. Intended User Profile:**

5.1. Only licensed clinicians are allowed to perform placement or removal of the ENT surgical sponge.

5.2. The clinicians shall select appropriate size of the ENT surgical sponge to

place into nasal cavity or ear canal. It is recommended to use one ENT surgical sponge on each side of the nose or ear. Place more than one ENT surgical sponges on each side of the nose or ear should be considered by clinicians.

5.3. The ENT surgical sponge shall be removal by clinicians within 72 hours after insertion.

5.4. The ENT surgical sponge to sale by, or on the order of, a physician.

#### **6. Contraindication:**

6.1. Absolute Contraindication:

Patients allergy to PVA (polyvinyl alcohol) materials.

6.2. Relative Contraindication:

Pregnant woman or the patients who have had any adverse reactions by using PVA related materials.

#### **7. Side Effects:**

7.1. An improper introduction and removal of the tamponade can lead to pain, bleeding and/or traumatization of the mucous membranes of the cavity inserted.

7.2. If the ENT surgical sponge is placed more than 72 hours, it might lead to local infection or systemic infection such as toxic shock syndrome.

7.3. No other side effects are known if applied properly.

#### **8. Precaution:**

8.1. Sterile, for single use only.

8.2. Sterility guaranteed unless unit package is opened or damaged.

8.3. Do not reuse, as it may result in compromised device performance and risk of cross contamination.

8.4. The ENT surgical sponge is not suitable to be moistened with an alcoholic solution as the structure of the sponge will be changed.

8.5. The ENT surgical sponge is not suggested to use on people who may be allergic to PVA sponge.

8.6. Please note that trimming the foam may cause risks, such as placement difficulty or displacement...etc.

8.7. Due to the unique production process of the tubing products, there is a chance that the following situations may occur:

As the product is stored in vacuum packaging, once the product is opened and moistened and has expanded, there is the possibility that the two ends of the silicon tube may not expand automatically. To avoid this condition, please carefully press and adjust the ends of the silicon tube before moistening the foam.

Please note that the above condition will not affect the original functionality of the product.

8.8. Dispose of used sponges according to local medical regulations.

### **9. Use Procedure:**

Please note that the PVA sponge will only expand properly by aqueous solutions moisten (do not use oil or fat-containing solutions/drugs to moisten).

Step 1. Prior to placing the ENT surgical sponge

- a. Select appropriate size of the ENT surgical sponge for placement into nasal cavity or ear canal.

Step 2. Removal of the ENT surgical sponge

- a. Prior to attempting retrieval of the ENT surgical sponge, the clinician may apply normal saline and wait for few, no more than 10 minutes until the ENT surgical sponge is thoroughly moistened.
- b. Using a tweezers to gently pull the ENT surgical sponge. If significant resistance is encountered, repeat step 2.a and pull again. Care must be taken not to damage surrounding mucosa if excessive adhesion occurs.
- c. Repeat step 2.b until the ENT surgical sponge is completely retrieved from the place inserted.

### **10. Warnings:**

- 10.1. Apply any medication onto the ENT surgical sponge for the purpose of hemostasis or anti-adhesion shall be directed or performed by a physician.
- 10.2. After applying the product more than three days, the evaluation by physician is suggested for further usage.
- 10.3. If any signs of irritation or infection to the nasal cavities or ear canal are observed, do not use the ENT surgical sponge and consult a physician immediately.
- 10.4. Rarely, displacement and ingestion of a nasal packing may occur. PVA surgical sponge is non-toxic and is considered safe to pass once it is completely in the alimentary tract. The patient should be monitored in the normal manner for the ingestion of any non-toxic material. Clinical confirmation that packing material has not been aspirated or is not in danger of being aspirated should be assessed. If aspiration should occur, prompt medical intervention is imperative with close patient monitoring and airway supporting pending the removal of the material.
- 10.5. In rare instances the physiochemical conditions associated with nasal or meatus or otitis chambers surgery, both with and without ENT surgical sponge, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include sudden fever, vomiting, diarrhea, dizziness, fainting (or near fainting when standing up), and/or a rash that looks like a sunburn. ENT

surgical sponge exhibit no antimicrobial properties, they are not bacteriostatic towards pre-existing infections, nor do they prevent the occurrence of new infections. In the case of pre-existing infections it is recommended that appropriate treatment be instituted.

**11. Preservation & Handling:**

- 11.1. The product should be preserved in a dry, well-ventilated place and in sealed conditions.
- 11.2. Do not use if package is opened or damaged and the identification tag with the Lot number should not be separated while in preservation.
- 11.3. Storage condition: 0°C~ 37°C.

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