Cenefom

Cenefom Corp.

3F., No.50-3, Keyan Rd., Hsinchu Science Park, Zhunan
Township, Miaoli County 35053, Taiwan

Tel: +886-37-580058

Fax: +886-37-580518

E-mail: service@cenefom.com www.cenefom.com

EC REP

MedNet EC-REP GmbH

Borkstrasse 10, 48163 Münster, Germany

Instructions for use

Ophthalmic Sponge

(for the series OP-ES, OP-EW, OP-CS, OP-LR and OP-LRD)

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English

1. Product Description:

The Ophthalmic Sponge is a sterile, single use device. The OP-ES series consisted of a spear is made of PVA (polyvinyl alcohol) material which can absorb aqueous fluids several times of its mass instantaneously upon contact with the fluids and a handle. The OP-EW, OP-CS, OP-LR and OP-LRD series are made of PVA (polyvinyl alcohol) material which can absorb aqueous fluids several times of its mass instantaneously upon contact with the fluids. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye).

Product Type	Applicable surgery or treatment
PVA Eye Spear or PVA Sponge Point	Absorbing excess fluids from the operative site.
Eye Wick	Drainage excess fluids and prevent them from entering the surgical site.
LASIK Ring	Absorbing and prevents the excess fluids and blood from encroaching into the surgical site.
LASIK Eye Drain	Absorbing excess fluids and protection of flap and hinge.
Corneal Light Shield	Keeping the cornea moist, protect the retina from potentially harmful light during ophthalmic procedures.

2. Intended Use:

The Ophthalmic sponge is indicated for use in ophthalmic or microscopic surgical procedures for absorbing excess fluids from the operative field.

3. Intended Therapeutic Use and Indications:

The Ophthalmic Sponge is mainly used for absorbing aqueous fluids in ophthalmic or microscopic surgical procedures.

4. Intended Patient Population:

Male and female older than 7 years.

5. Intended User Profile:

- 5.1. The ophthalmic sponge should be used on the order of a physician.
- 5.2. The ophthalmic sponge to sale by, or on the order of, a physician.

6. Contraindications

6.1. Absolute Contraindication:

Patient allergy to PVA (polyvinyl alcohol) materials.

6.2. Relative Contraindication:

Pregnant woman or the patients who have had any adverse reactions by using PVA related materials.

7. Side Effects:

No any side effects.

8. Precaution:

8.1. Sterile, for single use only.

8.2. Sterility guaranteed unless unit package is opened or damaged.

8.3. Do not reuse, as it may result in compromised device performance and risk of cross contamination.

8.4. The ophthalmic sponge is not suitable to be moistened with an alcoholic solution as the structure of the sponge will be changed.

8.5. The ophthalmic sponge is not suggested to use on people who may be allergic to PVA sponge.

8.6. Dispose of used sponges according to local medical regulations.

9. Use Procedure:

9.1. PVA Spear:

9.1.1. Opening *The Product Package*

9.1.1.1. The product package contains two foiled pouches sealed in an envelope. The envelope has a Tyvek top layer and a clear plastic bottom layer.

9.1.1.2. To open the package, separate the clear plastic layer from the Tyvek layer and gently pull open the envelope.

9.1.1.3. Empty the two foiled pouches onto a sterile field.
Caution: Maintain sterile technique when contacting the inner foiled pouches.

9.1.2. Opening *The Sterile Foiled Pouches*

9.1.2.1. Each foiled pouch contains five (5) PVA SPEAR.

9.1.2.2. While maintaining sterile technique, tear open the foiled pouch.

9.1.2.3. Remove spears and place in the sterile field.

Caution: In order to avoid inadvertent contact of the spear with aqueous fluids, keep the spears in a dry place within the sterile field.

9.1.3. *Using The PVA Spear*

- 9.1.3.1. By holding the plastic handle, gently touch the tip of the spear at the site where excess fluid is to be removed.
- 9.1.3.2. Hold still for approximately 3 seconds to allow for complete swollen of the PVA material.
Caution: Do not squeeze or press the PVA spear against the surface of the surgical sites.
- 9.1.3.3. Dispose of the used spear into a biological waste container.
- 9.1.3.4. When necessary, use additional spears to absorb fluid until all excess fluid is satisfactorily removed.

9.2. Eye Wick:

- 9.2.1. Moisten wick(s) with sterile fluid and insert around eye site (without pinching).
- 9.2.2. Irrigate normally.

9.3. LASIK Eye Drain & LASIK Ring:

- 9.3.1. Remove LASIK Eye Drain/LASIK Ring and place in the sterile field. Moisten LASIK Eye Drain/LASIK Ring if needed.
- 9.3.2. Use the LASIK Eye Drain/LASIK Ring at the surgical process site where excess fluid is to be removed.

9.4. Corneal Light Shield:

- 9.4.1. Remove Corneal Light Shield and place in the sterile field. Moisten Corneal Light Shield if needed.
- 9.4.2. Place the Corneal Light Shield onto the corneal to prevent light illumination and maintain cornea moist throughout surgery.

10. Warnings:

- 10.1. Do not attempt to use the sponge for cleaning of the surgical sites or removing debris, such as vitreous material that might be generated during the surgery.
- 10.2. The Ophthalmic sponge is not for use as an applicator for applying medication to the surgical sites. Also, the ophthalmic sponge shall not be used as a hemostasis tool for compressing the spear against the bleeding site.

11. Preservation & Handling:

- 11.1. The product should be preserved in a dry, well-ventilated place and in sealed conditions.
- 11.2. Do not use if package is opened or damaged and the identification tag with the Lot no. should not be separated while in preservation.
- 11.3. Storage condition: 0°C~ 37°C.

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