



# MEDICAL-GRADE PVA SPONGE

As a professional air-foaming Polyvinyl Alcohol (PVA) foam manufacturer, Cenefom provides high medical grade PVA Sponge for various applications and production purpose.



## BENEFITS



**Interconnected Pores**


- Interconnected pores structure results in extremely fast wicking
- Ultra- fine pore materials minimize adhesion
- Soft and comfortable when wet
- Lint and fiber-free
- Excellent tensile strength
- Well bio-compatibility
- Highly absorbent

## OEM/ ODM Service

We specialize at precise PVA foam supply according to customers' specific requirements.

### Physical Properties Spec.

Subjective	Apply Range
Pore size	200~1400 $\mu\text{m}$
Water absorption speed	min. 3 sec (17 mm)
Absorption capability	300~2000 %
Relative hardness	0.01°~85°
Thickness in wet	0.5~150 mm
Tensile strength in wet	min. 150 KPa
Formaldehyde residue	min. 2.5 ppm
Water Transmissibility	max. 96 %



# MEDICAL-GRADE PVA SPONGE

Cenefom General Medical PVA Sheet meets high medical grade standards. It is widely conducted in surgery fields like ophthalmic, ENT, wound dressing and medical tampon.

## PVA RAW SHEET



Model	Dimensions	Application
RS-01	270 x 180 x N* mm	Ophthalmic surgery
RS-02	270 x 180 x N* mm 490 x 250 x N* mm	Nasal packing
RS-03	330 x 220 x N* mm	NPWT dressing
RS-04	290 x 190 x N* mm	Incontinence tampon / other medical tampon

\*\* Height is available from 1 mm to 150 mm

Negative-pressure wound therapy (NPWT) is an advance therapeutic technique using a suction pump, tubing and a dressing to remove excess exudate and promote healing in acute or chronic wounds. Cenefom provides NPWT White Foam Dressing meets high standard requirements in wound care fields. With Cenefom's unique air-foaming PVA technology, we specialize at supplying high quality white foam dressing for NPWT usage.

## NPWT DRESSING

### Product Features

- Homogeneous distributed porosity with excellent drainage
- Strong porous structure to keep a smooth draining process
- Low-adhesion
- Moderate material softness with high standard cleanliness



Model	Dimensions	Packaging
RS-NT01	150 x 50 mm	20 pcs per bag in vacuum
RS-NT02	150 x 100 mm	10 pcs per bag in vacuum
RS-NT03	200 x 100 mm	10 pcs per bag in vacuum
RS-NT04	200 x 150 mm	10 pcs per bag in vacuum
RS-NT05	300 x 200 mm	5 pcs per bag in vacuum
RS-NT07	100 x 75 mm	20 pcs per bag in vacuum